




SimMan ALS



- Completely tetherless
- LinkBox and battery housed inside simulator
- Internal compressor (as in 3G)
- Limited only to battery life
- Unilateral chest rise
- Auto inflation (from SW):
 - Tongue Oedema
- Pneumothorax bladders

- Interchangeable pupils
- Manual moveable eyelids with three settings:
 - Open
 - Closed
 - 50%




- Bilateral carotid
- Bilateral femoral
- Right radial



Laerdal

helping save lives



- Anterior lung sounds
- Posterior lung sounds
- Heart sounds
- Left and Right Lung
 - Open / Medium / Closed
- Chest drain insertion
 - Bilateral

- Defibrillation uses ShockLink

**Purchased separately*

- Users **QCPR**
 - Rate
 - Depth
 - Hand position
 - Ventilation – rate / volume



Laerdal
helping save lives

- Cyanosis available
- Pre-ported IV cubital fossa

**No site for cannulation*

- Tibial IO
- Comes with male and female genitalia

**Dry procedure only*

Additional Ultrasound Features

- Comes Ultrasound prepared
- Will be able to be used with 3G scenario sets
- Will have EMS focused ultrasound scenario sets in the future